



The **GO2** is an adaptable O2 holder designed to safely and quickly transport medical oxygen for personal use in school bus, community transport, ambulance or non-emergency transportation vehicles.

- **GO2** prevents O2 cylinders from tipping, falling or becoming a projectile during a sudden stop or manoeuvre.
- **GO2** works with Series L Track or floor pocket for quick fleet-wide deployment and ideal vehicle and passenger placement.

SPECIFICATIONS:

O2 Cylinder Holder meets National School Transportation Specifications & Procedures 2015 IEP-IFSP Process, Guidelines, E.10


O2 Cylinder Holder meets Ambulance Manufacturers Division AMD Standard 003, Oxygen Tank Retention System

CEN Standard – BS EN 1789:2007, Medical Vehicles and their Equipment, 4.5.9 and 5.3



Ergonomic tightening and release knob

Includes two (2) anti-slip straps to hold O2 cylinder in place



Quick release that adjusts easily to accommodate O2 cylinders from sizes M2 to E in America as well as AZ to D in Europe

Made of high grade anodized aluminum for durability and corrosion resistance



Visual locking indicator to help ensure the GO2 is secured

Compact in size and weighs less than 2.72 Kg.



Mounts easily to existing **Series L Track or pockets¹** on vehicle floor and can be moved from one vehicle to another without the use of tools



The GO2 is also available in a Twin Cylinder Holder configuration.

This way, operators can accommodate multiple passengers that require O2 or that need a spare close at hand.

¹Excluding covered floor pockets



NEW!
TWIN HOLDER
AVAILABLE!



The **GO2** by **SURE-LOK** is an adaptable O2 holder designed to safely and quickly transport medical oxygen for personal use in school bus, para-transit, community transport, ambulance or non-emergency transportation vehicles.

- **GO2** prevents O2 cylinders from tipping, falling or becoming a projectile during a sudden stop or maneuver.
- **GO2** works with Series L Track for quick fleet-wide deployment and ideal vehicle and passenger placement.

SPECIFICATIONS:

O2 Cylinder Holder meets National School Transportation Specifications & Procedures 2015 IEP-IFSP Process, Guidelines, E.10

O2 Cylinder Holder meets Ambulance Manufacturers Division AMD Standard 003, Oxygen Tank Retention System

CEN Standard – BS EN 1789:2007, Medical Vehicles and their Equipment , 4.5.9 and 5.3

SURE-LOK 

O2 ON THE GO




The GO2 is available in its original single cylinder holder configuration as well as a **NEW TWIN cylinder holder.**

This way, operators can accommodate multiple passengers that require O₂ or that need a spare close at hand.



Ergonomic tightening and release knob

Includes two (2) anti-slip straps to hold O₂ cylinder in place




Quick release that adjusts easily to accommodate O₂ cylinders from M2 to E

Made of high grade anodized aluminum for durability and corrosion resistance



Visual locking indicator to help ensure the GO2 is secured

Compact in size and weighs less than 6 lbs.



Mounts easily to existing **Series L Track** or pockets¹ on vehicle floor and can be moved from one vehicle to another without the use of tools

¹Excluding covered floor pockets



THE **SURE-LOK GO2** OXYGEN (O₂) CYLINDER HOLDER IS DESIGNED TO SAFELY TRANSPORT MEDICAL OXYGEN FOR PERSONAL USE IN SCHOOL BUS, PARA-TRANSIT, COMMUNITY TRANSPORT, AMBULANCE OR NON-EMERGENCY TRANSPORTATION VEHICLES. IT HELPS TO PREVENT O₂ CYLINDERS FROM TIPPING, FALLING OR BECOMING A PROJECTILE DURING A SUDDEN STOP OR MANEUVER.

PRODUCT SPECIFICATIONS:

- Designed with tightening and release knob with visual lock indicator
- Will have two anti-slip straps to contain cylinder
- Shall be made of high grade anodized aluminum
- Stainless steel L Track fitting for easy installation and removal
- Compact design of 8" (203mm) x 6" (152mm) x 17" (432mm) and weighs less than 6 lbs. (2.72 Kg.)
- 1 Year Warranty

O₂ CYLINDER HOLDER FEATURES SHALL BE DESIGNED WITH:

- Ergonomic tightening and release knob for ease of use
- Visual lock indicator to help ensure the GO2 is secured
- Two anti-slip straps to hold O₂ cylinder in place while in transport
- Quick release for easy height adjustment to accommodate O₂ cylinders from size M2 to E (Americas) as well as AZ to D (Europe)
- High grade anodized aluminum for durability and resistance to corrosion
- Series L Track fitting mounts easily on vehicle floor and moves from one vehicle to another without the use of tools
- Compact size for ease of transportation and storage

TESTING SPECIFICATIONS:

- Dynamic crash test at 30mph/20g
- National School Transportation Specifications & Procedures 2015 IEP-IFSP Process, Guidelines, E.10
- Ambulance Manufacturer's Division AMD Standard 003, Oxygen Tank Retention System
- CEN Standard – BS EN 1789:2007, Medical Vehicles and their Equipment, 4.5.9 and 5.3

